Thermo Scientific Barnstead Hose Nipple Cartridges

Economical and simple way to purify water in small volumes up to 75L/hr



- Barnstead[™] Hose Nipple Cartridges may be installed in a series using simple tubing for greater capacity
- Color change indicator identifies when resin is exhausted in specific cartridges
- Designed for non-pressurized applications
- Cartridges include 0.38in. hose barb connection on each end except D8822, D8950 and D8951, which include tapered straight nipple for 0.38in. ID tubing
- Constructed of 100% virgin polypropylene cartridge with excellent quality resin

Applications: Softening, deionization, organic and chlorine removal **Alert:** Alcohol-containing samples cannot be used in cartridges containing the color indicator

System Warranty¥: One year parts, 90 days labor

Cartridge Warranty¥: 1 year parts

Cat. No.	Description	Flowrate	Length
D8900	Hose Nipple Cartridge Holder	Up to 75L/hr	1/2 Size: 25.9cm (10.2in.) 2/3 Size: 33.8cm (13.3in.) Full Size: 47.5cm (18.7in.)

Required Accessories

Cat. No.	Description	Resin Type(s)	Capacity	Color Indicator [†]
D50220	1/2 Size Mixed Bed	Ultrapure Ion Exchange	430 grains	No
D8822	1/2 Size Mixed Bed with Oxygen Removal	Ultrapure Mixed Bed with Strong Anion Exchange	280 grains/4.4g	No
D8950	2/3 Size High Capacity	Two bed Ion Exchange with Strong Cation Exchange	1100 grains	Yes
D8951	2/3 Size Mixed Bed and Organic Removal	Mixed Bed Ion Exchange with Activated Carbon	470 grains/ 2000 gal.	Yes
D8905	Full Size Cation Removal	Cation Exchange	3000 grains	Yes
D8901	Full Size High Capacity	Two Bed Ion Exchange	1650 grains	Yes
D8908	Full Size Macroreticular, Cation Removal	Macroreticular Cation Exchange	n/a	No
D8904	Full Size Organic Removal	Activated Carbon	5000 gal.	No
D8903	Full Size Oxygen Removal	Strong Anion Exchange	30g	No
D8921	Full Size Pretreatmentand Scale Eliminator	Mixed Bed and Strong Cation Exchange, Activated Charcoal	1250 grains/1000 gal.	Yes
D8911	Full Size Ultrapure	Ultrapure Ion Exchange	915 grains	No
D8902	Full Size Ultrapure	Ultrapure Mixed Bed Ion Exchange	915 grains	Yes
D8922	Full Size Ultrapure with Organic Removal	Mixed Bed Ion Exchange with Activated Carbon	730 grains/2000 gal.	Yes

†Alcohol-containing samples cannot be used in cartridges containing the color indicator.