

customizeable for your pretreatment or deionization needs

Harvey DI+ Cartridge System

- Turn-key system ships complete with a remote dispenser, Pura-Lite (50 KΩ·cm) indicator and a high capacity, two-bed ion exchange cartridge to remove calcium, magnesium, and other ions from water
- The Harvey DI+ system accepts a broad range of other cartridges to remove impurities and produces pure water at a flow rate of 4 L/min



| Product Specifications | | | | | |
|-------------------------|--|---------------------|---------------------------|----------------------------|------------------|
| | Overall dimensions W x H x D | Inlet connection | Feed water temperature | Max feed water pressure | Max flow rate |
| Single holder | 7 x 24 x 7 in (17.8 x 61 x 17.8 cm) | 1/2 in NPTF | 4-49°C (40-120°F) | 100 psig | 4 L/min |
| Double holder | 15 x 27 x 7 in (38.1 x 68.65 x 17.8 cm) | 1/2 in NPTF | 4-49°C (40-120°F) | 100 psig | 4 L/min |
| Half size holder | 7 x 15 x 7 in (17.8 x 38.1 x 17.8 cm) | 1/2 in NPTF | 4-49°C (40-120°F) | 100 psig | 2 L/min |
| Harvey DI+ | 7 x 26 x 7 in (17.8 x 63 x 17.8 cm) | 1/2 NPTF | 4-49°C (40-120°F) | 100 psig | 4 L/min |

Thermo Scientific Barnstead B-Pure (continued)



| Cartridges | |
|---------------------|------------------------------|
| Size | Dimensions Length x Diameter |
| Full Size Cartridge | 17.5 x 3.4 in |
| Half Size Cartridge | 10.95 x 3.4 in |

| System Options | Electrical | Model number |
|---|------------|--------------|
| Single holder | – | D4511 |
| Double holder with digital purity meter | 120 VAC | D4521 |
| | 240 VAC | D4522-33* |
| Double holder with Pura-Lite indicator (50 kΩ) | 120 VAC | D4524 |
| Double holder with Pura-Lite indicator (200 kΩ) | 120 VAC | D5831 |
| Double holder with Pura-Lite indicator (1 MΩ) | 120 VAC | D5833 |
| Half size holder | – | D4505 |
| Harvey DI+ | 120V | AY1273X4 |

*  